## Immunology: Mouse Neutralizing and Blocking Antibodies

## In vivo Functional Grade Cell-Specific Depletion Antibodies

Lot-to-lot reproducible monoclonal antibodies against mouse cell-specific proteins and cytokines are paramount for reproducible results. Depletion can occur naturally or be induced for treatment purposes. Therefore, it is a useful tool in therapeutic, immunotherapy and diagnostic research. It is common practice for researchers to use antibodies to remove a specific cell type to infer its function and mouse cell depletion *in vivo* experiments are used to establish the role of specific cell types in a variety of immunological processes. Many of these antibodies have been developed by academic and government labs all over the world and have been reported and referenced for use in *in vivo* studies to deplete a specific cell type.

## **Cell Types**

T Cells	Treg Cells	Monocytes
B Cells	γδ T Cells	Macrophages
CD8+ T Cells	NK Cells	Neutrophils
CD4+ T Cells	ILCs	Eosinophils

## In vivo GOLD™ & PLATINUM™ Quality Standards

*In vivo* GOLD™ and PLATINUM™ functional grade antibodies from Leinco Technologies are known for lot-to-lot reproducibility, low endotoxins and no detectable impurities.

Specifications	<i>in vivo</i> GOLD™ Mabs	<i>in vivo</i> PLATINUM™ Mabs			
Binding validation determined by WB, FC or ELISA	Yes	Yes			
Endotoxin level determined by the LAL method	≤ 1.0 EU/mg	≤ 0.5 EU/mg			
Antibody aggregation screening by analytical SEC	≥ 95% monomer	≥ 98% monomer			
Purity by SDS Page	≥ 95%	≥ 98%			
Formulation for <i>in vivo</i> use	No preservatives No stabilizers No carrier proteins Sterile PBS pH 7.2 - No K <sup>+</sup> or CA <sup>2+</sup> Concentration: > 5mg/mL	No preservatives No stabilizers No carrier proteins Sterile PBS pH 7.2 - No K <sup>+</sup> or CA <sup>2+</sup> Concentration: > 5mg/mL			
Murine pathogen screening	Not applicable, see PLATINUM functional grade antibodies	Pathogen tested (IMPACT1)			
Applications	In vivo functional studies and may be used for studies, as well as WB, FC, IF or IHC				
Product preparation	Functional grade preclinical antibodies are manufactured in an animal free facility using only <i>in vitro</i> protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.				
Storage and handling	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.  • 1 month, 2 to 8 °C under sterile conditions, as supplied.  • 12 months, -70 °C under sterile conditions.				





Scan the QR code to see our full catalog of Cell-Specific Depletion Antibodies

Cell Type	Target	Clone	GOLD Prod. No.	PLATINUM Prod. No.
T Cells	CD3ε	145-2C11	C1758	C2280
i Celis	Thy1.1 (CD90.1)	19E12	C3101	C3102
3	Thy1.2 (CD90.2)	30H12	C393	C430
	Thy (CD90)	T24/31	C392	C1789
3 X 5	TCR β	H57-597	T662	T772
	γ/δ	UC7-13D5	T352	T952
	CD4	YTS 191	C3220	C3210
<b>B Cells</b>	B220	RA3-6B2	C383	C2844
A. A.	CD19	1D3	C2117	C2849
	CD22	Cy34.1	C3358	C4358
				0 1000
CD8+T Cells	CD8a	2.43	C380	C2837
	CD8a	53-6.7	C375	C2848
	CD8a	YTS 169	C2442	C2850
	CD8β	53-5.8	C2832	C2836
	CD8(Lyt 2.1)	116 - 13.1	C3110	C3111
CD4+ T Cells				
	CD4	GK1.5	C1333	C2838
	CD4	YTS191	C3220	C3210
Treg Cells				
Y.	CD25 (IL-2Rα)	PC61	C1194	C2845
-	ODZS (IL-Zhu)	FC01	01194	02043
γδ T Cells				
4	Vγ2 TCR	UC3-10A6	T353	T704
NK Cells	NIIZ 4 4	DICLOS	N1100	Noon
	NK 1.1	PK136	N123	N268
	CD122 (IL-2Rβ)	TM-β1	C2325	C6325
ILCs	Thy1.2 (CD90.2)	30H12	C393	C430
Monocytes	Luco (0* 1)	DDC 005	G100	G153
	Ly6G (Gr-1) CSF1R	RB6-8C5 AFS98	C2169	C2268
	CSFIR	AF390	02109	02200
Macrophages				
	CSF1R (CD115)	AFS98	C2169	C2268
	CSFIR (CDII3)	AF596	C2169	U2200
Neutrophils	Ly6G	1A8	L280	1.005
	Ly6G (Gr-1)	RB6-8C5	G100	L305
	Lyou (GI-1)	ND0-0C3	G100	G153
Eosinophils				
(C)	IL-5	TRFK5	I-1060	I-1192
	IL-5	IHFNO	1-1060	1-1192
Plasmacytoid Dendritic Cells				
Denaritic Cells	CD317 (PDCA-1)	927	C791	C792
	/		-	
				410 Axminister Dr.



410 Axminister Dr. St. Louis, MO 63026 800.538.1145 leincoglobal@leinco.com www.leinco.com S01-2211-B